

FavorPrepTM **Tissue Total RNA Mini Kit**

User Manual

Cat. No.: FATRK 001 (50 Preps)
FATRK 001-1 (100 Preps)
FATRK 001-2 (300 Preps)

For Research Use Only

v.1101

Introduction

FavorPrep Tissue Total RNA Extraction Mini Kit is designed for extraction of total RNA from animal tissue and cultured cells. Some specially modified protocols are developed for other samples, such as bacteria and yeast. This method first lyses cells by using a chaotropic salt, then binds RNA to silica-based membranes, washes RNA with ethanol-contained wash buffer and then elutes purified RNA by RNase-free ddH₂O. It takes 30 min for an entire procedure, and the purified RNA is ready for RT-PCR, northern blotting, primer extension and cDNA library construction.

Sample amount and yield

Sample	Recommended amount of sample used		Yield (μg)
Animal cells (up to 5×10^6)	NIH/3T3	1×10^6 cells	10
	HeLa	1×10^6 cells	15
	COS-7	1×10^6 cells	30
	LMH	1×10^6 cells	12
Animal tissues (Mouse/rat) (up to 30 mg)	Embryo	10 mg	25
	Heart	10 mg	10
	Brain	10 mg	10
	Kidney	10 mg	30
	Liver	10 mg	50
	Spleen	10 mg	35
	Lung	10 mg	15
Thymus	10 mg	45	
Bacteria	<i>E. coli</i>	1×10^9 cells	60
	<i>B. subtilis</i>	1×10^9 cells	40
Yeast (up to 5×10^7)	<i>S. cerevisiae</i>	1×10^7 cells	25

Handling time: about 30 min

Kit Contents

Cat. No. / preps	FATRK001 (50 preps)	FATRK001-1 (100 preps)	FATRK001-2 (300 preps)
FARB Buffer	25 ml	45 ml	130 ml
Wash Buffer 1	30 ml	60 ml	175 ml
Wash Buffer 2 (concentrated)	15 ml *	35 ml **	50 ml *** x 2
RNase-free ddH ₂ O	6 ml	6 ml	8 ml X2
Filter Column	50 pcs	100 pcs	300 pcs
FARB Mini Column	50 pcs	100 pcs	300 pcs
Collection Tube	100 pcs	200 pcs	600 pcs
Micropestle	50 pcs	100 pcs	300 pcs
Elution Tube	50 pcs	100 pcs	300 pcs
User manual	1	1	1

* Add 60 ml ethanol (96-100 %) to Wash Buffer 2 when first open.


** Add 140 ml ethanol (96-100 %) to Wash Buffer 2 when first open.

*** Add 200 ml ethanol (96-100 %) to each Wash Buffer 2 when first open.

Important notes

1. Make sure everything is RNase-free when handling RNA.
2. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
3. Pipet a required volume of FARB Buffer to another RNase-free container and add 10 μ l β -mercaptoethanol (β -ME) per 1ml FARB Buffer before use.
4. Add required volume of RNase-free ethanol (96~100%) to Wash Buffer 2 as bottle indicated when first open.
5. All centrifuge steps are done at full speed (14,000 rpm or 10,000 x g) in a microcentrifuge.
6. Dilute RNase-free DNase 1 in reaction buffer (1M NaCl, 10mM MnCl₂, 20 mM Tris-HCl, pH 7.0 at 25°C) to final conc. = 0.5 U/ μ l.
7. The additional equipment, 20-G needle syringe, is needed for extraction of total RNA from tissue sample.

Brief Procedure



General Protocol: (For Animal Cells)

Please Read Important Notes Before Starting The Following Steps.

1. Pellet $1 \sim 5 \times 10^6$ cells by centrifuge at 300 x g for 5 min. Remove all the supernatant.
2. Add 350 μ l of FARB Buffer (β -ME added) to the cell pellet and vortex vigorously to lyse the cells. Incubate at room temperature for 5 min.

Note: In order to release all the RNA in the sample, it is required to disrupt the sample completely. Different samples require different methods (ex: disruptor equipment) to achieve complete disruption.

3. Place a Filter Column into a Collection Tube and transfer the sample mixture to Filter Column, centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 min.
4. Transfer the clarified supernatant from Collection Tube to a new micro-centrifuge tube (not provided) and adjust the volume of the clear lysate.
--Avoid pipetting any debris and pellet from Collection Tube.
5. Add 1 volume of 70% ethanol to the clear lysate and mix well by vortexing.
6. Briefly spin the tube to remove drops from the inside of the lid.
Place a FARB Mini Column into a Collection, transfer the ethanol added sample (including any precipitate) to FARB Mini Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min and discard the flow-through.
7. (Optional): To eliminate genomic DNA contamination, follow the steps from 7a. Otherwise, proceed to step 8 directly.
 - 7a. Add 250 μ l of Wash Buffer 1 to wash FARB Mini Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.
 - 7b. Add 100 μ l of RNase-free DNase 1 solution (0.5U/ μ l, not provided) to the membrane center of FARB Mini Column. Place the Column on the benchtop for 15 min.
 - 7c. Add 250 μ l of Wash Buffer 1 to wash FARB Mini Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.
 - 7d. After DNase 1 treatment, proceed to step 9.
8. Add 500 μ l of Wash Buffer 1 to wash FARB Mini Column. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.

9. Wash FARB Mini Column twice with 750 μ l of Wash Buffer 2 by centrifuge at full speed (14,000 rpm or 10,000 x g) for 1 min then discard the flow-through.
--Make sure that ethanol has been added into Wash Buffer 2 when first open.
10. Centrifuge at full speed (14,000 rpm or 10,000 x g) for an additional 3 min to dry the column.
--**Important Step!** This step will avoid the residual liquid to inhibit subsequent enzymatic reaction.
11. Place FARB Mini Column to Elution Tube.
12. Add 50 μ l of RNase-free ddH₂O to the membrane center of FARB Mini Column. Stand FARB Mini Column for 1 min.
--**Important Step!** For effective elution, make sure that RNase-free ddH₂O is dispensed on the membrane center and is absorbed completely.
13. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 min to elute RNA.
14. Store RNA at -70C.

Special Protocol: (For Animal Tissue)

Additional equipment: a 20-G needle syringe

1. (For Fresh sample): Cut up to 30 mg of tissue sample. Grind the tissue sample completely in liquid nitrogen. Transfer the powder to a new microcentrifuge tube(not provided). Or you can place tissue sample into a microcentrifuge tube and use provided micropestle to grind the tissue sample few times and break it into small pieces.
(For Frozen sample) Weight up to 30 mg tissue sample and grind the tissue sample in liquid nitrogen then transfer the powder to a new microcentrifuge tube (not provided).
2. Add 350 μ l of FARB Buffer (β -ME added) to the sample and shear this tissue sample by passing lysate through a 20-G needle syringe 10 times. Incubate at room temperature for 5 min.
--Grind the sample a few times to make it break more completely.

Note: In order to release all the RNA in the sample, it is required to disrupt the sample completely. Different samples require different methods (ex: disruptor equipment) to achieve complete disruption.

3. Follow the Animal Cells Protocol starting from step 3.

Special Protocol: (For Bacteria)

1. Transfer 1 ml well-grown bacterial culture (or up to 1×10^9 cells) to a microcentrifuge tube (not provided).
2. Descend the bacterial cells by centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 min and discard the supernatant completely.
3. Resuspend the cell pellet in 100 μ l of RNase-free lysozyme reaction solution (20mg/ml lysozyme; 20mM Tris-HCl, pH 8.0; 2mM EDTA; 1.2% Triton) (not provided).
4. Incubate at 37°C for 10 min.
5. Add 350 μ l of FARB Buffer (β -ME added) to the sample and mix well by vortex. Incubate at room temperature for 5 min.

Note: In order to release all the RNA in the sample, it is required to disrupt the sample completely. Different samples require different methods (ex: disruptor equipment) to achieve complete disruption.

6. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 min to spin down insoluble material and transfer the supernatant to a microcentrifuge tube (not provided) and adjust the volume of the clear lysate.
-Avoid pipetting any debris and pellet in the Collection Tube.
7. Follow the Animal Cells Protocol starting from step 5.

Special Protocol: (For Yeast)

1. Transfer 3 ml of log-phase ($OD_{600}=10$) yeast culture to a microcentrifuge tube (not provided).
2. Descend the yeast cells by centrifuge at 7,500 rpm (5,000 x g) for 10 min and discard the supernatant completely.
3. Resuspend the cell pellet in 600 μ l sorbitol buffer (1 M sorbitol; 100 mM EDTA; 0.1% β -ME) (not provided). Add 200 U zymolase or lyticase and incubate at 30 °C for 30 min.
--Prepare sorbitol buffer just before use.
4. Centrifuge at 7,500 rpm (5,000xg) for 5 min to pellet the spheroplasts. Discard the supernatant.
5. Add 350 μ l of FARB Buffer (β -ME added) to the pellet and vortex vigorously to lyse the spheroplasts. Incubate at room temperature for 5 min.

Note: In order to release all the RNA in the sample, it is required to disrupt the sample completely. Different samples require different methods (ex: disruptor equipment) to achieve complete disruption.

6. Centrifuge at full speed (14,000 rpm or 10,000 x g) for 2 min to spin down insoluble materials and transfer the supernatant to a microcentrifuge tube (not provided).
7. Follow the Animal Cells Protocol starting from step 5.